



SRP14 Antibody, HRP conjugated

Product Code	CSB-PA022673LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37108
Immunogen	Recombinant Human Signal recognition particle 14 kDa protein (2-107AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal recognition particle 14 kDa protein (SRP14) (18 kDa Alu RNA-binding protein), SRP14
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SRP14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.