



TCERG1 Antibody, HRP conjugated

Product Code	CSB-PA023291LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14776
Immunogen	Recombinant Human Transcription elongation regulator 1 protein (782-896AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription elongation regulator 1 (TATA box-binding protein-associated factor 2S) (Transcription factor CA150), TCERG1, CA150 TAF2S
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TCERG1