



TMEM88 Antibody, Biotin conjugated

Product Code	CSB-PA023892LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PEY1
Immunogen	Recombinant Human Transmembrane protein 88 protein (109-159AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 88, TMEM88, TMEM88A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM88
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.