



TMF1 Antibody, Biotin conjugated

Product Code	CSB-PA023904LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P82094
Immunogen	Recombinant Human TATA element modulatory factor protein (878-986AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TATA element modulatory factor (TMF) (Androgen receptor coactivator 160 kDa protein) (Androgen receptor-associated protein of 160 kDa), TMF1, ARA160
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMF1