



TMF1 Antibody, Biotin conjugated

Product Code	CSB-PA023904LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P82094
Immunogen	Recombinant Human TATA element modulatory factor protein (878-986AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TATA element modulatory factor (TMF) (Androgen receptor coactivator 160 kDa protein) (Androgen receptor-associated protein of 160 kDa), TMF1, ARA160
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.