



TSPAN12 Antibody, Biotin conjugated

Product Code	CSB-PA025149LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95859
Immunogen	Recombinant Human Tetraspanin-12 protein (111-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-12 (Tspan-12) (Tetraspan NET-2) (Transmembrane 4 superfamily member 12), TSPAN12, NET2 TM4SF12
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TSPAN12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.