



LTN1 Antibody, Biotin conjugated

Product Code	CSB-PA026646LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94822
Immunogen	Recombinant Human E3 ubiquitin-protein ligase listerin protein (462-591AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase listerin (EC 2.3.2.27) (RING finger protein 160) (RING-type E3 ubiquitin transferase listerin) (Zinc finger protein 294), LTN1, C21orf10 C21orf98 KIAA0714 RNF160 ZNF294
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LTN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.