



KLRC3 Antibody, FITC conjugated

Product Code	CSB-PA23099C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07444
Immunogen	Recombinant Human NKG2-E type II integral membrane protein (94-240AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NKG2-E type II integral membrane protein (NK cell receptor E) (NKG2-E-activating NK receptor), KLRC3, NKG2E
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KLRC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.