



KLRC3 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA23099D0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q07444 |
| Immunogen | Recombinant Human NKG2-E type II integral membrane protein (94-240AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | NKG2-E type II integral membrane protein (NK cell receptor E) (NKG2-E-activating NK receptor), KLRC3, NKG2E |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | KLRC3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |