



PKD1 Antibody, HRP conjugated

Product Code	CSB-PA23149B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98161
Immunogen	Recombinant Human Polycystin-1 protein (615-753AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycystin-1 (PC1) (Autosomal dominant polycystic kidney disease 1 protein), PKD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PKD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.