



TTLL3 Antibody, Biotin conjugated

Product Code	CSB-PA23329D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4R7
Immunogen	Recombinant Human Tubulin monoglycylase TTLL3 protein (470-641AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin monoglycylase TTLL3 (EC 6.3.2.-) (HOTTL) (Tubulin--tyrosine ligase-like protein 3), TTLL3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TTLL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.