



# KIF23 Antibody, HRP conjugated

Product Code	CSB-PA23569B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02241
Immunogen	Recombinant Human Kinesin-like protein KIF23 protein (865-960AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF23 (Kinesin-like protein 5) (Mitotic kinesin-like protein 1), KIF23, KNSL5 MKLP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KIF23