



# KIF23 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA23569C0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q02241
<b>Immunogen</b>	Recombinant Human Kinesin-like protein KIF23 protein (865-960AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Kinesin-like protein KIF23 (Kinesin-like protein 5) (Mitotic kinesin-like protein 1), KIF23, KN5L5, MKLP1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	KIF23
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.