



KIF23 Antibody, Biotin conjugated

Product Code	CSB-PA23569D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02241
Immunogen	Recombinant Human Kinesin-like protein KIF23 protein (865-960AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF23 (Kinesin-like protein 5) (Mitotic kinesin-like protein 1), KIF23, KNLS5 MKLP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KIF23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.