



PKIG Antibody, Biotin conjugated

Product Code	CSB-PA23619D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2B9
Immunogen	Recombinant Human cAMP-dependent protein kinase inhibitor gamma protein (1-76AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	cAMP-dependent protein kinase inhibitor gamma (PKI-gamma), PKIG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PKIG