




FPGS Antibody

Product Code	CSB-PA23719A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05932
Immunogen	Recombinant Human Folylpolylglutamate synthase, mitochondrial protein (443-546AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Folylpolylglutamate synthase, mitochondrial (EC 6.3.2.17) (Folylpolyl-gamma-glutamate synthetase) (FPGS) (Tetrahydrofolylpolylglutamate synthase) (Tetrahydrofolate synthase), FPGS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FPGS

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, MCF-7 whole cell lysate

All lanes: FPGS antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 65, 61, 60, 62 kDa

Observed band size: 65 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.