



FPGS Antibody, HRP conjugated

Product Code	CSB-PA23719B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05932
Immunogen	Recombinant Human Folylpolyglutamate synthase, mitochondrial protein (443-546AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Folylpolyglutamate synthase, mitochondrial (EC 6.3.2.17) (Folylpoly-gamma-glutamate synthetase) (FPGS) (Tetrahydrofolypolyglutamate synthase) (Tetrahydrofolate synthase), FPGS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FPGS