



HLA-F Antibody, FITC conjugated

Product Code	CSB-PA335552LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30511
Immunogen	Recombinant Human HLA class I histocompatibility antigen, alpha chain F protein (194-303AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class I histocompatibility antigen, alpha chain F (CDA12) (HLA F antigen) (Leukocyte antigen F) (MHC class I antigen F), HLA-F, HLA-5.4 HLAF
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-F
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.