



MPP3 Antibody, HRP conjugated

Product Code	CSB-PA615679LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13368
Immunogen	Recombinant Human MAGUK p55 subfamily member 3 protein (301-432AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MAGUK p55 subfamily member 3 (Discs large homolog 3) (Protein MPP3), MPP3, DLG3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MPP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.