



IQGAP2 Antibody, FITC conjugated

Product Code	CSB-PA615701LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13576
Immunogen	Recombinant Human Ras GTPase-activating-like protein IQGAP2 protein (1245-1395AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras GTPase-activating-like protein IQGAP2, IQGAP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IQGAP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.