



SGCG Antibody, HRP conjugated

Product Code	CSB-PA618771LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13326
Immunogen	Recombinant Human Gamma-sarcoglycan protein (108-221AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-sarcoglycan (Gamma-SG) (35 kDa dystrophin-associated glycoprotein) (35DAG), SGCG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SGCG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.