



TRIM26 Antibody, HRP conjugated

Product Code	CSB-PA619638LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12899
Immunogen	Recombinant Human Tripartite motif-containing protein 26 protein (151-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 26 (EC 2.3.2.27) (Acid finger protein) (AFP) (RING finger protein 95) (Zinc finger protein 173), TRIM26, RNF95 ZNF173
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRIM26
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.