



KIR2DS5 Antibody, Biotin conjugated

Product Code	CSB-PA619902LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14953
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 2DS5 protein (162-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell immunoglobulin-like receptor 2DS5 (CD158 antigen-like family member G) (MHC class I NK cell receptor) (Natural killer-associated transcript 9) (NKAT-9) (CD antigen CD158g), KIR2DS5, CD158G NKAT9
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KIR2DS5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.