



# QPRT Antibody, Biotin conjugated

Product Code	CSB-PA621868LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15274
Immunogen	Recombinant Human Nicotinate-nucleotide pyrophosphorylase [carboxylating] protein (113-254AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nicotinate-nucleotide pyrophosphorylase [carboxylating] (EC 2.4.2.19) (Quinolinate phosphoribosyltransferase [decarboxylating]) (QAPRTase) (QPRTase), QPRT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	QPRT