



DAZAP2 Antibody, FITC conjugated

Product Code	CSB-PA622986LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15038
Immunogen	Recombinant Human DAZ-associated protein 2 protein (1-168AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DAZ-associated protein 2 (Deleted in azoospermia-associated protein 2), DAZAP2, KIAA0058
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DAZAP2