



PTPN14 Antibody, HRP conjugated

Product Code	CSB-PA623010LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15678
Immunogen	Recombinant Human Tyrosine-protein phosphatase non-receptor type 14 protein (220-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein phosphatase non-receptor type 14 (EC 3.1.3.48) (Protein-tyrosine phosphatase pez), PTPN14, PEZ PTPD2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTPN14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.