



PLCL1 Antibody, FITC conjugated

Product Code	CSB-PA624015LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15111
Immunogen	Recombinant Human Inactive phospholipase C-like protein 1 protein (1-108AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive phospholipase C-like protein 1 (PLC-L1) (Phospholipase C-deleted in lung carcinoma) (Phospholipase C-related but catalytically inactive protein) (PRIP), PLCL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLCL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.