



MPEG1 Antibody, FITC conjugated

Product Code	CSB-PA642654LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2M385
Immunogen	Recombinant Human Macrophage-expressed gene 1 protein (537-633AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Macrophage-expressed gene 1 protein (Macrophage gene 1 protein) (Mpg-1), MPEG1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MPEG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.