



MPEG1 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA642654LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q2M385 |
| Immunogen | Recombinant Human Macrophage-expressed gene 1 protein (537-633AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Macrophage-expressed gene 1 protein (Macrophage gene 1 protein) (Mpg-1), MPEG1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | MPEG1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |