



FREM1 Antibody, Biotin conjugated

Product Code	CSB-PA686013LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5H8C1
Immunogen	Recombinant Human FRAS1-related extracellular matrix protein 1 protein (1957-2082AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FRAS1-related extracellular matrix protein 1 (Protein QBRICK), FREM1, C9orf143 C9orf145 C9orf154
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	FREM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.