



HYI Antibody, HRP conjugated

Product Code	CSB-PA719395LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T013
Immunogen	Recombinant Human Putative hydroxypyruvate isomerase protein (1-174AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative hydroxypyruvate isomerase (EC 5.3.1.22) (Endothelial cell apoptosis protein E-CE1), HYI
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HYI
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.