



IBA57 Antibody, Biotin conjugated

Product Code	CSB-PA722566LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T440
Immunogen	Recombinant Human Putative transferase CAF17, mitochondrial protein (286-356AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative transferase CAF17, mitochondrial (EC 2.1.-.-) (Iron-sulfur cluster assembly factor homolog), IBA57, C1orf69
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	IBA57