



IBA57 Antibody, Biotin conjugated

Product Code	CSB-PA722566LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T440
Immunogen	Recombinant Human Putative transferase CAF17, mitochondrial protein (286-356AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative transferase CAF17, mitochondrial (EC 2.1.-.-) (Iron-sulfur cluster assembly factor homolog), IBA57, C1orf69
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	IBA57
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.