



SLC6A19 Antibody, HRP conjugated

Product Code	CSB-PA731656LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q695T7
Immunogen	Recombinant Human Sodium-dependent neutral amino acid transporter B(0)AT1 protein (326-413AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent neutral amino acid transporter B(0)AT1 (Solute carrier family 6 member 19) (System B(0) neutral amino acid transporter AT1), SLC6A19, B0AT1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC6A19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.