



SLC6A19 Antibody, Biotin conjugated

Product Code	CSB-PA731656LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q695T7
Immunogen	Recombinant Human Sodium-dependent neutral amino acid transporter B(0)AT1 protein (326-413AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent neutral amino acid transporter B(0)AT1 (Solute carrier family 6 member 19) (System B(0) neutral amino acid transporter AT1), SLC6A19, B0AT1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC6A19