



VAR2 Antibody, FITC conjugated

Product Code	CSB-PA735951LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5ST30
Immunogen	Recombinant Human Valine--tRNA ligase, mitochondrial protein (830-981AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Valine--tRNA ligase, mitochondrial (EC 6.1.1.9) (Valyl-tRNA synthetase) (ValRS) (Valyl-tRNA synthetase-like), VAR2, KIAA1885 VAR2L VARSL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VAR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.