



RIPPLY2 Antibody, Biotin conjugated

Product Code	CSB-PA735979LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TAB7
Immunogen	Recombinant Human Protein ripply2 protein (1-76AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ripply2, RIPPLY2, C6orf159
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RIPPLY2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.