



THEMIS2 Antibody, HRP conjugated

Product Code	CSB-PA735984LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TEJ8
Immunogen	Recombinant Human Protein THEMIS2 protein (62-177AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein THEMIS2 (Induced by contact to basement membrane 1 protein) (Protein ICB-1) (Thymocyte-expressed molecule involved in selection protein 2), THEMIS2, C1orf38 ICB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	THEMIS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.