



# THEMIS2 Antibody, Biotin conjugated

Product Code	CSB-PA735984LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TEJ8
Immunogen	Recombinant Human Protein THEMIS2 protein (62-177AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein THEMIS2 (Induced by contact to basement membrane 1 protein) (Protein ICB-1) (Thymocyte-expressed molecule involved in selection protein 2), THEMIS2, C1orf38 ICB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	THEMIS2