



SLC27A1 Antibody, Biotin conjugated

Product Code	CSB-PA740768LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PCB7
Immunogen	Recombinant Human Long-chain fatty acid transport protein 1 protein (453-570AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Long-chain fatty acid transport protein 1 (FATP-1) (Fatty acid transport protein 1) (EC 6.2.1.-) (Solute carrier family 27 member 1), SLC27A1, ACSVL5 FATP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC27A1