



SLC27A1 Antibody, Biotin conjugated

Product Code	CSB-PA740768LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PCB7
Immunogen	Recombinant Human Long-chain fatty acid transport protein 1 protein (453-570AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Long-chain fatty acid transport protein 1 (FATP-1) (Fatty acid transport protein 1) (EC 6.2.1.-) (Solute carrier family 27 member 1), SLC27A1, ACSVL5 FATP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC27A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.