



MPZL3 Antibody, Biotin conjugated

Product Code	CSB-PA740916LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UWV2
Immunogen	Recombinant Human Myelin protein zero-like protein 3 protein (32-158AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin protein zero-like protein 3, MPZL3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MPZL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.