



GSG1L Antibody, FITC conjugated

Product Code	CSB-PA740924LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXU4
Immunogen	Recombinant Human Germ cell-specific gene 1-like protein (84-176AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Germ cell-specific gene 1-like protein (GSG1-like protein), GSG1L
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GSG1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.