



PARS2 Antibody, Biotin conjugated

Product Code	CSB-PA745334LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L3T8
Immunogen	Recombinant Human Probable proline--tRNA ligase, mitochondrial protein (126-258AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable proline--tRNA ligase, mitochondrial (EC 6.1.1.15) (Prolyl-tRNA synthetase) (ProRS), PARS2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PARS2