



MEX3B Antibody, HRP conjugated

Product Code	CSB-PA751152LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZN04
Immunogen	Recombinant Human RNA-binding protein MEX3B protein (253-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein MEX3B (RING finger and KH domain-containing protein 3) (RING finger protein 195), MEX3B, KIAA2009 RKHD3 RNF195
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MEX3B