



MEX3B Antibody, HRP conjugated

Product Code	CSB-PA751152LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZN04
Immunogen	Recombinant Human RNA-binding protein MEX3B protein (253-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein MEX3B (RING finger and KH domain-containing protein 3) (RING finger protein 195), MEX3B, KIAA2009 RKHD3 RNF195
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MEX3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.