



DPY19L2 Antibody, FITC conjugated

Product Code	CSB-PA757514LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NUT2
Immunogen	Recombinant Human Probable C-mannosyltransferase DPY19L2 protein (599-749AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable C-mannosyltransferase DPY19L2 (EC 2.4.1.-) (Dpy-19-like protein 2) (Protein dpy-19 homolog 2), DPY19L2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DPY19L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.