



UBR3 Antibody, Biotin conjugated

Product Code	CSB-PA757959LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZT12
Immunogen	Recombinant Human E3 ubiquitin-protein ligase UBR3 protein (960-1113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase UBR3 (EC 2.3.2.27) (N-recognin-3) (RING-type E3 ubiquitin transferase UBR3) (Ubiquitin-protein ligase E3-alpha-3) (Ubiquitin-protein ligase E3-alpha-III) (Zinc finger protein 650), UBR3, KIAA2024 ZNF650
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.