




TFAP2E Antibody

Product Code	CSB-PA761516LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6VUC0
Immunogen	Recombinant Human Transcription factor AP-2-epsilon protein (66-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor AP-2-epsilon (AP2-epsilon) (Activating enhancer-binding protein 2-epsilon), TFAP2E
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TFAP2E

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate

All lanes: TFAP2E antibody at 3.63µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 47 kDa

Observed band size: 47 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.