



HOGA1 Antibody, HRP conjugated

Product Code	CSB-PA771464LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86XE5
Immunogen	Recombinant Human 4-hydroxy-2-oxoglutarate aldolase, mitochondrial protein (1-164AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	4-hydroxy-2-oxoglutarate aldolase, mitochondrial (EC 4.1.3.16) (Dihydrodipicolinate synthase-like) (DHDPS-like protein) (Probable 2-keto-4-hydroxyglutarate aldolase) (Probable KHG-aldolase) (Protein 569272), HOGA1, C10orf65 DHDP5L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HOGA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.